
FJ-4000 MASSAGE CHAIR

OPERATION INSTRUCTIONS

SAVE THESE INSTRUCTIONS

Thanks for purchasing this product. Please read this manual instruction carefully in order to know how to use this product correctly before operation. Please pay more attention to some important safety warning and keep this manual book properly for future use.

Note: Our Company reserves the right to revise the design and description of this product without any further notice. Please refer to the actual color of this product.

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Safety and maintenance

1. Important Safety Warning.

- Don't allow children to touch moveable part of this product.
- Please use well-grounding power supply suited to this product.
- Please pull out the plug after use or before cleaning to avoid injury or Damage to the product.
- Please operate this product according to this manual instruction.
- Don't use accessories which are not recommended.
- Don't use this product outdoors.
- Please read this manual instruction carefully before operation.
- Any other usage that is not listed in this manual instruction is forbidden.
- 20 minutes of usage each time is recommended.
- Please don't use it if the leather is broken or damaged.
- Please don't use it if the cover or leather is damaged.
- Please don't use it if the radiator hole is covered or blocked.
- Please don't drop anything into this product.
- Please don't fall asleep while using this product.
- Please don't use it if drunk or feel unwell.
- Please don't use it within one hour after having meal.
- Please don't make the massage function too strong for avoiding Injury.

2. Environment for usage.

- Please don't use it under high temperature and moist environment such as bathroom.
- Please don't use it immediately while the environmental temperature changes sharply.
- Please don't use it under heavy-dusty or caustic environment.
- Please don't use it where there is not enough space or not good ventilation.

3. People who are not suitable to use this product.

- People who are suffering from osteoporosis are not suitable to use this product.
- People who are suffering from heart disease or carrying electronic medical equipment such as heart pacemaker are not suitable to use this product.
- People who are suffering from fever are not suitable to use this product.
- Pregnant women or women that are in menses period are not suitable to use this product.
- People who got injured or skin disease are not suitable to use this product.
- Children under 14-year old or people mentally unbalanced without being Supervised are not allowed to use this machine.
- People who are told by doctor to have rest or people who feel unwell are not suitable to use this product.
- People whose bodies are wet are forbidden to use this product.
- The appliance has a heated surface. persons insensitive to Heat must be careful when using the appliance.

4. Safety.

- Check the voltage if it is suitable to the specification of this product.
- Don't use wet hand to pull the plug.
- Don't make water go inside this product to avoid electric shock or cause damage to this product.

Safety and maintenance

- Don't damage the wires or change the circuit of this product.
- Don't use wet cloth to clean the electric parts such as switch and plug.
- Be away from this product under power cut state to avoid injury if the power comeback suddenly.
- Stop using this product while it works abnormal, and please consult local agent immediately.
- Stop using this product if you feel unwell, and please consult health care doctor.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.
- If the supply cord is damaged, it must be replaced by the manufacturer, its Service agent or similarly qualified persons in order to avoid a hazard.

5. Maintenance.

- This product should be maintained by certain appointed agents, users are prohibited to disassemble or maintain by themselves.
- Please don't forget to cut the power after usage.
- Don't use this product if the socket looses.
- If this product will be left unused for long, please curl up the wires and reserve this product in dry and dustless environment.
- Don't reserve this product under high temperature or near fire and avoid longtime direct sunshine.
- If the detachable cord / supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a Hazard.
- Please clean this product by dry cloth. Don't use thinner, benzene or alcohol.
- The mechanical components of this product are specially designed and made, it is not necessary to be specially maintained.
- Don't use edged res to thrust this product.
- Don't roll or pull this product above uneven ground, it should be lifted before moving.
- Please use it intermittently; don't make this product work continuously for long.


6. Solution of ordinary malfunctions.

- It is nature movement voice if the motor sounds while in usage.
- If the controller can't work normally, please check the plug and socket if they are connected strongly and make sure the switch is turned on.
- If the rated working time is over, the product will automatically activate the power-off button; If this product works continuously for long, the temperature protector will make the product automatically power-off, it should be used again after half an hour's rest.

Use safety warning


⚠ Warning

- When the footrest lifting, people should not apply pressure on the footrest to avoid the massage chair falling forward to damage the product or cause personal injuries

⚠ Warning


- Before sitting on the massage chair, Please make sure the massage hand cannot place in the seat position, so as avoiding the massage hand to damage the product or cause personal buttocks injuries.

Use safety warning


⚠ Warning

- Do not strive to pull out the control cable or power cord during standby or power-on state to avoid product damaged or electrocution accidents. Advise you unplug the controller cable and power cord and put away after you use the chair.


⚠ Warning


- When turning on the zero-gravity lying function, be sure to notice whether there are children or pets under the massage table under the bottom of the calf and the backrest to stay or play to avoid product damage and personal injury.

Name and function of components

- 1.Head
- 2.Shoulder air bag components
- 3.Manual controller pocket
- 4.Arm air bag components
- 5.Seat cushion
- 6.Calf rest
- 7.Foot rest
- 8.Calf rest end sealing plate
- 9.Back cushion
- 10.Armrest
- 11.Armrest decoration
- 12.Side cover


- 13.Back cover
- 14.Remote controller
- 15.PCB box cover
- 16.Power cord
- 17.Controller cord socket
- 18.Switch
- 19.Fuse box
- 20.Power cord socket

Function

- This product is designed with a set of four- wheel driven mechanical massage hands that can move up and down while stretching back and forth.
- Designed with five emulation massage methods of shiatsu, kneading, flapping, knocking, simultaneous kneading and flapping
- Set with six kinds of characterized auto massage function: comfort, ease, relax, neck& shoulder massage, back& waist massage, fast improving experience massage.
- Set with 1 memory function
- Manually select the type of upper body massage with the five massage point options which include fixed position, partial position, full back, neck& shoulder, back& waist. With three level of speeds for adjustable. In state of fixed position, the massage ball can move up and down to adjust the massage point accurately. In the state of tapping, shiatsu and knocking, the width between the two kneading balls is adjustable with three levels: wide, medium and narrow.
- Air pressure massage function with three modes: full body air pressure, arms air pressure, leg air pressure. with 3 intensity levels adjustable.
- Leg designed with roller massage function. At the bottom of foot, two pairs of papillae on toe root rotating massage; on the arch and heel mastoid skin scraping massage.
- The back heating function uses carbon fiber as an infrared heating source.
- The footrest can be extended to suit different leg lengths.
- set the calf and the back of a chair lift regulating function alone. The calf will be auto separate the mechanism when it touch foreign body.
- zero space, reach the maximum space saving, and design with three auto-relining mode.
- L shape back track fitting your body curves , can realize the zero gravity when back down
- The chair will be reset when switch off. The massage hands will be reset as well.


Method of usage

Note: The specific displays refer to the actual object of this product!

Method of usage


1.Switch power supply on

2.Start massage

- 1).Press power button, start massage.
- 2).Automatic shoulder position and massage acupressure point detection, waiting for the completion of testing.
- 3).If there is difference with the massage balls and the actual shoulder position, you can press the up and down button to adjust the ball to a suitable position. Shoulder position has 10 grades adjustable. It will run automatic massage if there is no adjustment on shoulder position within 10 seconds.

Button	Description	Display
	Move the massage balls upward to a suitable position	
	Move the massage balls downward to a suitable position	

4).Choose massage function to massage

· Memory function massage: press memory it will revert to the previous memory massage mode after “Di”. (Totally three memory mode, M1, M2, M3), otherwise " DI,DI", indicating that this mode hasn't been stored before.

Function	Button	Description	Display
Memory Massage		Memory Massage includes three modes	Memory

Method of usage

• Characterized auto massage

Function	Button	Description	Display
Memory Massage		Ache relieve: Deeply massage on acupressure points, effectively alleviate pain	Auto F1
		Fatigue recovery: Depth massage, awaken body energy	Auto F2
		Comfortable: Loosen the muscles.	Auto F3
		Stretch: The calf-rest and back-rest lifting together, combine with air pressure massage to stretch the calf ,which relieve fatigue and revitalizing	Auto F4
		Fast: Experience all massage functions quickly	Auto F5
		Buttock: Focused on buttock massage	Auto F6
		Neck & shoulder: Focused on neck and shoulder massage	Auto F7
		Back & waist: Focused on back and waist massage	Auto F8

Note:  The massage includes four kinds of auto modes: ache relieve, fatigue recovery, comfort and stretch.

 The massage includes two kinds of auto modes: shoulder and back & waist

• Upper body manual massage

Function	Button	Description	Display
Technique		Shiatsu : Speed And Width	
		Knead & flap: Speed	
		Flap: Speed And Width	
		Knead: Speed	
		Knock: Speed And Width	
Speed		Three Grades Speed Adjustable	
Width		Three Grades Width Adjustable	

Method of usage

Function	Button	Description	Display
Back Stretch		Fixed position (fixed spot) massage	 H0
		Small area (partial area) massage back and forth	 H1
		Overall massage back and forth	 H2
		Upper back massage back and forth	 H3
		Lower back massage back and forth	 H4
		Massage point adjust (in fixed spot, partial modes you can adjust the massage position)	Massage hand move up
			Massage hand move down

• Air pressure function massage

Function	Button	Description	Display
Air pressure mode		Three air pressure modes (full body, arms and legs) for lower body	Arm  Feet 
Intensity		With three grades adjustable	Arm  Feet 

• Other manual function

Function	Button	Description	Display
Roller		Sole roller: sole roller switch, with 3 speed levels available	
Heating		Far infrared back heating, has warm feeling when it start to work 3 minutes later	

• Other manual function

Function	Button	Description
Calf rise		Calf rest lifting button, press this button, calf rest lifting slowly, stop when release it
Calf drop		Calf rest decline button, press this button, calf rest decline slowly, stop when release it
Backrest rise		Backrest lifting button, press this button, backrest lifting slowly, stop when release it
Backrest drop		Backrest decline button, press this button, backrest decline slowly, stop when release it
Auto lay down mode		With three modes circulation function, when press the button, turn into auto lay down mode.

Method of usage


3. Turn off power, stop massage

- 1) Press power button during massage progress, all massage function will be stop immediately, the backrest and footrest will return to original position automatically, meanwhile you can set up memory function. After rated working time, all massage functions will be finished, meanwhile you can set up memory function, but backrest and footrest won't return to original position.
- 2) Memory store function: After rated working time or switch off, the hand controller displays "Memory" and glitter (The display will prompt 1, 2, 3, if you want to store a memory mode, just stay in this mode for two seconds until the hand controller alarming "Di".) , remind the users whether need to store the massage mode before turn off.
- 3) Cut off the power of the whole machine. Figure (Switch power off)


4. Calf installation

- The install hole of rail shelf (chart 1).
- Place the isolation sleeve on both sides of the calf hook. (Outside) (chart2).


Method of usage

- Place the flat washer into the pin shaft. (chart3).
- Calf hooks install hole and rail shelf install hole should be aligned, insert the pin with a flat washer(both sides are required to be inserted). The pins are aligned with the screw holes for easy installation. (chart4).


- Use the Allen wrenches lock two hexagon socket head screws. (chart5).

5. Head cushions and backrest pad's dismounting and use manual


- Foldable pillow pad. The users can choose the number of layer to find the best intensity of back massage.
- Separate the head pillow to fall down to the back to reduce the intensity of back massage.
- Backrest cushion and backrest connected by zipper (1), head cushion and backrest pad connected by thread gluing (2)

Method of usage

6.Floor Protection


- Put the heavy massage chairs on the wood floor for a longtime may damage the floor, so please place a carpet or other object to avoid this damage.

Caution: Do not move the chairs with its wheels on wood floor, rough ground or in a narrow space. Two people are required to lift the chair. During the lifting period, take the calf sealing board and backrest cover grooves as the focal point(Please do not lift armrest Directly.)

7.Method of movement

Press the remote button to make the calf to be lifting to the highest level, then turn off the power switch, unplug the power wire, remote controller wire. Hold a groove in the back plastic cover to tilt the backrest backward to a certain degree (with the gravity center resting on the wheels), push the chair forward or backward with your hands and finally resume the chair to the normal position in a slow and gentle manner.

Note:During movement you must turn off the power switch first and then unplug the power line and hand controller connecting line.

Product specification

Model: FJ-4000

Name: Massage Chair

Rated voltage: 110-120V~ 60Hz ☐
220-240V~ 50Hz/60Hz ☐

Rated power input:150W

Rated time: 20Min

Safety structure: Class I

Material introduction: PVC, PA, steel parts and electric & electron parts

MATERIAL	PU	LEATHER	CLOTH	WOOD

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